




Nicovape Q user guide

read before use

Welcome to your Nicovape Q




*sold separately


Getting started

1. Remove the coloured safety cap and mouthpiece plug from your Q cartridge

 Please recycle



2. Insert the Q cartridge into the top of your Nicovape Q device. It will secure itself automatically via the magnets



Getting started

3. For the best experience, ensure your device is fully charged by using the USB-C charger provided

The light ring around the charge port will indicate your device's current charge level



No charge



Running low



Normal charge



Fully charged

A flashing white light ring means your Q Cartridge may be empty, or there may be a connection error. Try removing it and re-inserting, or replace with a new cartridge if it appears to be empty.

4. Once your device is fully charged, please refer to the reverse side of this booklet for important guidance on how to start using your new device

How to get the best experience


The Nicovape Q is simple and familiar to use, and activates only when you draw through the mouthpiece.

- 1.** To start, place the mouthpiece between your lips as shown on the right, about halfway along the mouthpiece.
- 2.** Draw firmly on the mouthpiece for 2 to 3 seconds - Aim to draw enough air through the device to fill your mouth.
- 3.** Immediately inhale deeply and exhale slowly.


Wait at least 20 seconds before your next puff. This will enable you to see if you've had enough to satisfy your cravings and will allow the device to reset for the next puff.

It's not unusual to cough when you first start, but this will ease with continued use.

Just stick with it!



Ensure half of mouthpiece is
between lips



Ensure device is pointing
downwards

Nicovape Customer Care

Mon - Fri 9am - 5pm AEST



Customer Care hotline

1300 642 682



Live online chat

nicovapeq.com